

# STERIS V-PRO 60



## ABOUT

---

The **Steris V-PRO 60** low temperature sterilizer is designed to provide reliable sterilization for heat- and moisture-sensitive medical devices. Using vaporized hydrogen peroxide (VHP) technology, the V-PRO 60 ensures effective sterilization while maintaining compatibility with a wide range of materials. This compact system offers three pre-programmed cycles—Lumen, Non-Lumen, and Flexible—to address diverse sterilization needs. Its automated processes, safety features, and rapid cycle times make it an efficient choice for modern healthcare facilities.

## FEATURES

---

**Low Temperature Sterilization:** Operates at low pressure and temperature, ideal for heat- and moisture-sensitive devices.

**VHP Technology:** Utilizes VAPROX HC Sterilant to produce vaporized hydrogen peroxide for thorough sterilization.

**Three Pre-Programmed Cycles:** Lumen Cycle: for devices with diffusion-restricted spaces; Flexible Cycle: for flexible endoscopes or bronchoscopes, Non-Lumen Cycle: for non-lumened rigid, semi-rigid, and flexible devices.

**Automated Processes:** Fully automated sterilization cycle with condition, sterilize, and aeration phases for precision and reliability.

**Compact Design:** Available in single-door counter or cart-mounted configurations to fit various facility layouts.



## DIMENSIONS AND WEIGHT

---

**Sterilizer Dimensions:**

31" W x 31" L x 28" H (78.7 x 78.7 x 71.1 cm)

**Chamber Dimensions:**

13" W x 28" L x 10" H (33.0 x 71.1 x 25.4 cm)

**Chamber Volume:**

2.1 cubic feet (60 L)

**Weight:**

120 Vac units: 255 lb (116 kg)

230 Vac units: 259 lb (117 kg)

# SPECIFICATIONS



## STERILIZATION CYCLES

The V-PRO 60 Low Temperature Sterilization System offers three pre-programmed cycles, each tailored for specific sterilization needs:

**Lumen Cycle:** Approximately 60 minutes, designed for devices with lumens or diffusion-restricted spaces.

**Non-Lumen Cycle:** Approximately 28 minutes, for sterilizing non-lumened instruments.

**Flexible Cycle:** Approximately 38 minutes, intended for one flexible endoscope or similar devices.

## CONDITION PHASES

Each cycle progresses through three phases: Condition, Sterilize, and Aeration.

**Condition Phase:** During this phase, the sterilizer removes air and moisture from the chamber. The reservoir fills, and a vacuum pulse is applied. The load is tested for moisture content; if acceptable, the process advances. If moisture levels are too high, the pulse repeats until conditions are met.

**Sterilize Phase:** This phase involves four pulses, each consisting of:

- Pulling a vacuum to the setpoint.
- Introducing vaporized VAPROX HC Sterilant into the chamber.
- Holding the vapor for a programmed time.
- Introducing filtered air to reach a setpoint, followed by another hold time.
- Finally, a deep vacuum is applied to complete the pulse.

**Aeration Phase:** In the final phase, a vacuum removes remaining sterilant vapors, reducing their concentration within the chamber. Once completed, chamber pressure returns to normal, and the door unlocks, signaling the end of the cycle.

## DISPLAY CONTROLS AND PRINTING

**Display Control Panel** is centrally located on the front of the sterilization unit for user convenience. This 5.7-inch color touch screen, with a resolution of 640 x 480 pixels, provides clear information and easy navigation. The interface is intuitive, featuring color-coded screens to streamline operation and enhance usability.

**Printer** is positioned on the front right side of the unit. It features an alphanumeric impact design and produces a clear, permanent record of each sterilization cycle. The printer uses 2.25-inch (5.7 cm) wide, 24-character tape, and includes a paper take-up for organized documentation. This system ensures accurate and accessible sterilization records for compliance and review.