



The CUSA EXcel+ system builds upon the market-leading features of the CUSA EXcel System with the addition of several performance enhancing elements. The CUSA EXcel+ system supports a broad selection of handpieces and surgical tips, allowing the surgeon to select the most appropriate combination based on the surgical procedure.

The user-friendly console, color coded tubing, and tubing paths enable easy system set-up. The CUSA EXcel+ system offers a simple user interface, the widest variety of handpieces and tip combinations, and advanced features that make it the most powerful ultrasonic tissue ablation system available.

## Compact

- Compact, easy to maneuver unit takes up minimal OR space

## Simple

- Simplified handpiece and manifold tubing connections.
- Easy, step-by-step console setup cards with color-coded function buttons.
- Steam sterilization reduces turnaround time and is in compliance with industry standards.

## Quiet

- Pumps run at lower values to reduce noise in the OR.
- Enters a quiet “resting” state when not in use.



# Specifications

**Dimensions**

Height: 132cm  
Width: 55.8cm  
Depth: 71.12cm  
Weight: 95.3kg

**Power**

Voltage 220-240V  
Power 110-120 VA

**Liquid**

Water Unmineralized water for cooling system  
Vacuum On board pump

**Cooling System**

The device is equipped with a closed water cooling system

**Environmental**

Functioning temperature 10 to 40°C  
Hygrometry 15 to 80% of humidity, without condensation

