




The **Axia STL4** is a mobile imaging table which is intended to support a 500 pound patient. The Axia STL4 is sold new and is intended to be used with mobile x-ray systems for minimally invasive procedures. The STL4 is made in the USA and has 4 table movements. This table is available with 2 different style table tops; standard and cervical. The cervical top is 22 inches wide with a fully adjustable headrest.

Features

- Easy c-arm access with diving board style table top.
- Large patient capacity up to 500 pounds.
- 2-inch thick mattress pad.
- Facial cutout for prone procedures.
- Easy patient transfer with low height.
- Available with a standard or cervical top.



Standard Top



Cervical Top

Specifications

Tabletop Dimensions **Length:** 84" (214.4 cm)
 Width: 24" (61cm)
 Thickness: 1.62" (4.1 cm)

Radiolucent Dimensions **Length:** 84" (142.2 cm)
 Width: 24" (61cm)

Table Dimensions **Length:** 86" (218.4 cm)
 Width: 26" (66 cm)

Weight **Table Weight:** 590lbs (268kg)
 Patient Capacity: 500lbs (227kg)

Table Movements **Height:** 30" to 40" (76.2 to 101.6 cm)
 Lateral Tilt: (+ / -) 17°
 Trendelenburg: 13°
 Longitudinal Travel: 10" (+/- 5) (25.4 cm (+/- 12.7 cm))

General **Mattress Thickness:** 2" (5 cm)
 Tabletop / Mattress Attenuation: 1.6mm Al.(+/- .2mm tolerance)
 Tabletop Material: Carbon fiber skin/foam core with integral facial cut-out
 Casters: 5" casters with full swivel and wheel lock
 Steel Base: I-shaped with epoxy finish
 Hand Control: Activates all motions

Power 120V/60HZ or 230V/50HZ
 Protect product from moisture
 Allowable Storage & Operating Humidity: 0-92%
 Allowable Storage Temperature: 32°F to 104°F (0°C to 40°C)
 Allowable Operating Temperature: up to 104°F (40°C)

