

Surgical Table

STERIS AMSCO 2080L / 2080RC

ABOUT

This Steris surgical table is available in two different models: the **Steris Amsco 2080L** and **2080RC**. Designed for a wide range of surgical procedures, this table is equipped with a removable headrest, a user-friendly hand pendant control unit (on the RC model), and a set of electrically conductive tabletop cushions for enhanced safety and comfort. Powered by a reliable 120V/60Hz electrical system, it ensures seamless operation and efficiency in the operating room. With a weight capacity of 300 lb, this table provides stable and secure support for patients.



FEATURES

Pedestal includes tabletop lift cylinder and support column with bearing-mounted saddle frame fully enclosed by stainless-steel, telescoping shrouds.

Each shroud is a one-piece construction to ensure that foreign matter enters the elevating mechanism.

300 lb weight capacity,

Side rails for accessories.

Adjustable headrest.

Movements include back up/down, leg up/down, Trendelenberg/reverse Trendelenberg, kidney elevator control, and head section control.

The Tabletop is a 20-inch wide tabletop divided into four hinged sections: head, back, seat, and foot.

Each section is cast aluminum (back, seat and foot sections are covered by stainless-steel veneer).

SPECIFICATIONS

DIMENSIONS

Height: 27 - 45 in (68.6 - 113.4 cm)

Width (back/seat plate): 20 in (50.8 cm)

Width (with side rails): 22 in (55.9 cm)

Length (w/o headrest): 71.25 in (181 cm)

Length (with headrest): 84 in (211.7 cm)

Weight: 690 lb (313 kg)

WEIGHT CAPACITY

Maximum Load (Lift): 300 lb (136 kg)

Maximum Load (Articulation): 300 lb (136 kg)

TABLE MOVEMENTS

Height: 27 - 45 in (68.6 - 113.4 cm)

Lateral Tilt (left/right): 10°

Trendelenburg: ± 20°

Reverse Trendelenburg: ± 20°

Kidney Bridge: 4 - 0.5 in (11.3 cm)

Back Up/Down: +90° / -10°

Leg Up/Down: +0° / -105°

Flex: -15° back/ -20° seat

Reflex: +10° back/ +15° seat