

STILLE SONESTA 6202

ABOUT

The **Stille Sonesta 6202** urology procedure chair offers optimal positioning for your facility. It features three movements and provides a wide range of optional accessories for various procedures. The Sonesta 6202 simplifies positioning, particularly for elderly or physically challenged patients. This lightweight chair can be easily moved, supports up to 375 pounds, and is easy to clean and disinfect after each procedure.

FEATURES

Three motors that control the height, seat, and back adjustments.

Height adjusts from 16 to 39.4 inches.

Lightweight design with locking mechanisms for secure positioning.

Supports patients weighing up to 375 pounds.



SPECIFICATIONS

DIMENSIONS

Height (with seat cushion): 16 in - 37 in (40 cm - 95 cm)

Width: 25 in (62 cm)

Length (with extension cushion): 75 in (190 cm)

Length (without extension cushion): 56 in (143 cm)

Weight: 276 lb (125 kg)

GENERAL

Lift Capacity: 375 lb (170 kg)

Wheel Diameter (Locking): 4 in (10 cm)

MOVEMENTS

Motorized Height, seat, and back adjustments

PROCEDURES

Urodynamics

General

Examinations

Cystoscopy

Ultrasound

Biofeedback

Vasectomy

OB/GYN Procedures

ELECTRICAL

120/230v 24v DC Output Meets EMC (Electromagnetic Compatibility) for IEC 60601 All Products carry the CE mark.