

TUTTNAUER 2340M



ABOUT

The Tuttnauer 2340M Autoclave Sterilizer is engineered to meet all sterilization needs at an affordable price while ensuring quality, safety, and reliability. This manual autoclave is packed with features for both safety and user convenience. It operates quietly and includes heat insulation to maintain a comfortable environment during operation.

FEATURES

Automatic shutoff at the end of both sterilization and drying cycles streamlining operations and conserving energy, while enhancing workflow efficiency.

Durable electro-polished chamber and door for extended use.

High-quality materials ensure long-lasting performance and maintain sterilization effectiveness over time.

Double safety lock to prevent the door from opening under pressure, prioritizing user safety, minimizing the risk of accidental exposure.

Built-in protection safeguards the autoclave, extending the life of the heating components.

Front-located drain valve for easy and quick water reservoir drainage

Dual safety thermostat to prevent overheating, this added layer of protection enhances safety during operation and promotes equipment longevity.

Complies with recognized standards, ensuring reliable and high-quality sterilization results.



SPECIFICATIONS

DIMENSIONS

Height: 14.4 in (365 mm)

Width: 20 in (510 mm)

Depth: 21.5in (545 mm)

Weight: 85 lbs (38.5 kg)

CHAMBER DIMENSIONS

Diameter: 9 in (230 mm)

Depth: 18 in (470 mm))

TRAYS

Number of trays: 3

Height:0.8 in (20 mm)

Depth: 16.3 in (415 mm)

Width: 6.7 in (170 mm)

Tuttnauer 2340M

SPECIFICATIONS



CASSETTE

Standard Cassettes Capacity: 2 full / 2 half.

ELECTRICAL

Voltage: 120V.

Frequency: 50 / 60 Hz.

Current: 12 A.

Power: 1400 W.

STANDARD UNWRAPPED CYCLE TIME

Cold: 27 Min.

Hot: 13 Min.

