

# SKYTRON 3600B ULTRASLIDE



## ABOUT

---

The Skytron 3600B UltraSlide surgical table offers a 23-inch Top Slide feature, delivering full table functionality and superior patient positioning for optimal C-arm access. With a robust 1,000 lb. lift capacity and 800 lb. articulation, it supports a wide range of procedures. Its removable back and leg sections, along with the Auto Beach Chair position, provide versatile patient positioning options. The UltraSlide 3600B is equipped with a durable, back-lit, one-touch pendant control for easy operation and enhanced visibility in low-light settings.



## FEATURES

---

Supports a lift capacity of up to 1,000 lbs. and an articulation capacity of 800 lbs., accommodating diverse surgical procedures.

A 23-inch top slide feature provides full table functionality, enabling superior patient positioning and maximizing C-arm access.

Includes a removable back and leg section with an Auto Beach Chair position for enhanced flexibility in patient setup.

A 57-inch imaging window allows comprehensive full-body imaging without requiring special accessories.

Features a durable, back-lit, one-touch pendant control for intuitive operation and clear visibility in low-light environments.

Minimum table height of 24 inches supports ergonomic and comfortable surgical access.

Battery power operation and a convenient return-to-level feature enhance reliability and ease of use..

## SPECIFICATIONS

### DIMENSIONS

---

**Height:** 24-42.5 in (61-108 cm)

**Width (Without side rails):** 19.5 in (50 cm)

**Width (With side rails):** 22 in (56 cm)

**Length (Without headrest):** 69 in (175 cm)

**Length (With headrest):** 82 in (208 cm)

**Weight:** 800 lbs (365 kg)

### MAXIMUM LOADS

---

**Lift:** 1,000 lb.

**Articulation:** 800 lb.

# SPECIFICATIONS



## MOVEMENTS

---

Height: 24-42.5 in (61-108 cm)

Lateral tilt (left/right): 30°

Trendelenburg: 30°

Reverse Trendelenburg: 30°

Power Kidney Bridge: 5.75 in lift (14.5 cm)

Back (up / down): 90°/40°

Leg down: 15°/100°

Return to level (anatomically correct return to horizontal position from articulated positions)

Top (slide/head): 9.5 inches 24 cm

Top (slide/foot): 13.5 inches 34 cm

Flex / Reflex: 40°/90°

Base braking (lock / unlock)

Beach Chair