

OAKWORKS CFPM100

ABOUT

The **Oakworks CFPM100** one-movement pain management table features a single tower lift system with a 500-pound weight capacity. Lowering to a height of 22-inches, the CFPM100 offers safe and secure wheelchair transfers. It's a cost-effective 22" x 84" carbon fiber table with 56" of metal-free imaging space. This CFPM table from Oakworks has an integrated hand rail for easy maneuverability.

FEATURES

- Minimum of 22-inches height range
- 500-pound weight capacity
- Carbon fiber table top
- 56-inches of metal-free imaging space factor.



SPECIFICATIONS

DIMENSIONS

Height: 22 – 40 in (56 – 102 cm)

Rectangular Table Top (WxL): 24 x 80 in (61 x 203 cm)

Integrated Headrest (WxL): 22 x 84 in (56 x 214 cm)

Weight: 343 lb (156 kg)

GENERAL

Patient Weight Capacity: 500 lb (227 kg)

Padding: 2 in (5 cm)

Casters: 4 in (10 cm)

Foot Control: Height motion

POWER

Voltage: 120V 60Hz or 230V/50Hz