

AXIA UROMAX 4

ABOUT

The **Axia UroMax 4** urology C-arm table offers an excellent solution for real-time urological X-ray procedures. The cantilevered, carbon fiber, radiolucent top areas allow full fluoroscopic visualization and unobstructed C-arm positioning. UroMax 4 has highly versatile movements, options, and accessories. The carbon fiber leg extension accessory can quickly turn your urology table into a procedure table for pain management and other C-arm-guided procedures.

FEATURES

- Carbon fiber with an integral facial cutout section
- Leg extension easy to remove and very lightweight
- Diving board style for easy C-arm access
- Large patient capacity (500-lb pound load)
- 2" patient high-quality comfort pad
- Light weight – easy to physically move the table
- 24" wide top eliminates the need for arm boards
- Easy patient transfer to table
- Low height for patient boarding
- Stirrups are an optional accessory



SPECIFICATIONS

DIMENSIONS

Table Length: 55 in (139 cm)

Table Width: 24 in (60.9 cm)

Table Height: 30 – 40 in (76.2 – 101.6 cm)

Weight Capacity: 500 lbs (226.7 kg)

Table Weight: 740 lbs (335.7 kg)

H-shaped Steel Base: 51 x 26 in (129.5 – 66 cm)

Mattress: 2 in thick (5 cm)

Radiolucent Area: 27 x 24 in (68.5 x 60.9 cm)

Table Top: 55 x 24 x 1.62 in (139 x 60 x 4 cm)

TILTS

Lateral Tilt: $\pm 17^\circ$

Trendelenburg: $\pm 13^\circ$

Longitudinal Travel: 10 in (± 5 in)

POWER

AC Voltage: 120V~

Frequency: 60 Hz

AC Current: 6.5A

Electric Shock Protection: Class I

Applied Part Rating: Type B

Mode of Operation: 10% Duty Cycle; 1 min on, 9 min off

ENVIRONMENTAL

Storage Humidity: 0- 92%

Operating Humidity: 0- 92%

Storage Temperature: 32°F – 104°F (0°C – 40°C)

Operating Temperature: Up to 104°F (40°C)