


AXIA HDVIEW D-SERIES LARYNGOSCOPE

ABOUT

The **Axia HDView D-Series** video laryngoscopes adopt state-of-the-art image technology and possess the advantage of a 2.0-megapixel full-view camera with a high-resolution monitor. **HDView D-Series** features Disposable blades, a rechargeable Lithium-ion battery for extended use and, a unique anti-fog capability is also provided. The overall design is compact and ergonomic. A high-resolution view of the glottis will enable clinicians first attempt success while minimizing any chance of a complication during the intubation process.

FEATURES

- 3.5 In Display
- Anti Fog Lens
- 2.0 Mega Pixel Camera
- Lightweight
- Anti-Shatter Protection Screen
- Disposable Blades



SPECIFICATIONS

DIMENSIONS

Height: 8.1 in (cm)
Width: 6 in (cm)
Depth: 2.5 in (cm)
Weight: 0.8 lb (350 g)

DISPLAY

Size: 3.5 In Full View
Monitor: LCD
Resolution: 640 x 480 RGB
Aspect Ratio: 4:3
Video Refresh Rate: 30 FPS
Illumination: LED

DATA OUTPUT

Easy to establish and store files

IMAGE OUTPUT

Connect to external monitor and convenient for training and presentation

LARYNGOSCOPE BLADE CAMERA

Mega Pixels: 2
Field Angle: 66°
Illuminance: = 800 LUX

WORKING CONDITIONS

Temperature: 5 ~?40
Humidity: 30 ~ 85
AP: 700 hPa to 1060 hpa

BATTERY

Type: Rechargeable lithium battery
Long Working Time: = 200 minutes
Voltage: 3.7V
Capacity: 3200 mAh
Battery Life Cycle: More than 300 charges
Charging Time: Less than 8 hours

TRANSPORT / STORAGE CONDITIONS

Temperature: -10 ~ ?45
Humidity: = 93
AP: 500 hPa to 1060 hPa

HANDLE

Comfortable ergonomic short handle design