




The **Axia SL1** was designed for high patient throughput for surgery centers, pain management, and other multidisciplinary medical facilities. The SL1 is a mobile imaging table which is intended to support a 500 pound patient. This Axia imaging table is intended to be used with mobile x-ray for minimally invasive procedures. This table is available with 2 different style table tops; standard and cervical. The cervical top is 22 inches wide with a fully adjustable headrest.

Features

- Carbon fiber, diving board style for easy C-arm access
- Adjustable height 28" - 38"
- Large patient capacity – high-end actuator motor supports 500-pound loads
- An 84" length can accommodate even the tallest of patients
- A 24" top eliminates the need for arm boards
- Two-inch thick, high quality patient comfort pad
- Facial cut-out for prone procedures
- Available with a standard or cervical top



Standard Top



Cervical Top

Specifications

Dimensions

Height (*Adjustable*): 28 - 38 In (71.1 - 96.5 cm)

Width: 24 In (61 cm)

Length: 86 In (218.4 cm)

Weight: 375 lbs (170 kg)

Tabletop

Height: 2 In (5 cm)

Width: 24 In (61 cm)

Length: 84 In (213.4 cm)

Material: Carbon fiber skin / foam core with integral facial cut-out

Radiolucent Area

Width: 24 In (61 cm)

Length: 56 In (142.2 cm)

General

Hand Control: Activates all motions

Base: I-Shaped with Epoxy Finish

Castors: 5 In (12.7 cm) with full swivel & wheel lock

Patient Straps: 2

