

MIDMARK RITTER 222 BARRIER FREE



ABOUT

The **Midmark Ritter 222 Barrier Free** powered exam table is designed to lower to just 18 inches from the floor, reducing the risk of injury from lifting patients. The Midmark Ritter 222 has a manual gas spring assist back, compared to the Midmark Ritter 223 which has a fully powered adjustable back. The Ritter 222 examination table makes your patients feel more stable and relaxed while waiting for their consultation as their feet rest on the floor.

FEATURES

- **Lowers to 18 inches from the floor**, providing easy access and transfer of patients even if they're in wheelchairs.
- **Maximum height of 37 inches**, elevating the patients to a wide array of height options will enable health professionals to conduct examinations with optimal posture and positioning.
- **Precise and smooth patient positioning with minimal effort**
- **Built-in pillows** makes the Ritter 222 a great platform for long and comfortable examinations.
- **Built-in stirrups that lock in four lateral positions** supports patients legs and prevent fatigue from awkward positioning.
- **Foot extension shelves** expand the table's length enabling the patient to lay down or rest their feet.
- **Foot control** provides quick and easy height control without any physical effort.



SPECIFICATIONS



DIMENSIONS

Height: 18" to 37" (45.7 to 94 cm)

Width: 30.5" (77.5 cm)

Length (Without foot extension): 59.5" (151.1 cm)

Length (With foot extension): 76" (193 cm)

Weight: 325 lbs (147 kg)

Patient Weight Capacity: 400 lbs (181.4 kg)

RANGE OF MOTION

Height: 18" to 37" (45.7 to 94 cm)

Back: 0° to 80° ± 5° (Gas spring assist)

Pelvic Tilt: Max: 7°; Min: 1.5°

PAPER ROLL HOLDER

Holds 21" x 3.5" paper rolls (53.3 cm x 8.9 cm)

ELECTRICAL

115 VAC, 50/60 Hz, 11 Amps (with heater & receptacle –016)

230 VAC, 50/60 Hz, 2.75 Amps (basic unit –017, excludes receptacle)

SAFETY STANDARDS

UL 60601-1

CAN/CSA No. 601.1-M90