

HAEMONETICS CELL SAVER ELITE+



ABOUT

The **Haemonetics Cell Saver Elite+** is an autotransfusion device designed to collect and recycle a patient's blood during and after surgery, allowing it to be processed, washed, and returned to the patient. Autotransfusion systems like the Cell Saver Elite+ are especially valuable in surgeries with significant blood loss. Soma Tech Intl provides refurbished Haemonetics Cell Saver Elite+ devices, expertly restored in Soma's ISO:13485-certified facility. Note that Soma Tech Intl operates independently and is not affiliated with Haemonetics.

FEATURES

- Reduces the need for allogeneic transfusions with autologous transfusion.
- Packs and washes red blood cells with over a 90% recovery rate.
- Touch-screen color monitor.
- Customizable settings.

RED BLOOD CELL QUALITY

- Lab testing of the 225 mL bowl with default settings.
- HCT: 60%
- RBC Recovery: 94%
- WBC Removal: 24.7%
- Free Hemoglobin Washout: 98.8%
- Albumin Washout 97.7%
- Heparin Washout 99.6%
- Potassium Washout: 96.4%
- Heparin Washout: 99.8%
- Fat Washout: 99.6%



DIMENSIONS

- Height: 48 in (121.9 cm)
- Width: 21 in (53.3 cm)
- Depth: 26.5 in (67.3 cm)
- Weight: 95 lb (43 kg)