

BIODEX ULTRA PRO

MODEL 056

ABOUT

The **Biodex Ultra Pro** ultrasound table is a specialized table for all general ultrasound procedures. With an ergonomic design catering to patients and sonographers, this table enhances the overall efficiency of achieving high-quality images. The open frame design facilitates closer proximity for sonographers, which is particularly beneficial for tasks like carotid imaging or long-leg arterial scanning.

FEATURES

Motorized height

Manual fowler positioning

Foot Support

Optional trendelenburg motion

Optional folding siderails

Optional foot control

Optional central floor locking system secures table position and is accessible from either side of the table



*Image based on certain options

SPECIFICATIONS

DIMENSIONS

Height (Adjustable): 23 in to 39 in (58.4 to 99 cm)

Width: 30 in (76.2 cm)

Length: 70 in (177.8 cm)

Weight: 335 lb (151.9 kg)

Patient Weight Capacity: 350 lb (158.7 kg)

5" (12.7 cm) Individual locking swivel casters

MATTRESS

Torso section 2" (5 cm) thick, center and leg sections 3" (8 cm) thick;

Foam rubber pad covered with seamless, durable Naugahyde

Upholstery Color: Graphite

MOTIONS

Height: 23 in to 39 in (58.4 to 99 cm)

Optional Trendelenburg: 0° to ±15°

Fowler Back (Manual): 0° to 80°

CONTROLS

Hand Control: Activates height and optional Trendelenburg motions

Manual fowler positioning

Foot Control (optional)

BIODEX ULTRA PRO

MODEL 058

ABOUT

The **Biodex Ultra Pro** ultrasound table is a specialized table for all general ultrasound procedures. With an ergonomic design catering to patients and sonographers, this table enhances the overall efficiency of achieving high-quality images. The open frame design facilitates closer proximity for sonographers, which is particularly beneficial for tasks like carotid imaging or long-leg arterial scanning. With the optional motorized Fowler back, this table offers infinite adjustability up to 80°.

FEATURES

The central floor locking system secures table position and is accessible from either side of the table

Motorized height and trendelenburg

Foot support

Manual Fowler positioning (Optionally motorized)

Optional folding side rails

Optional foot control



SPECIFICATIONS

DIMENSIONS

Height (Adjustable): 23 in to 39 in (58.4 to 99 cm)

Width: 30 in (76.2 cm)

Length: 70 in (177.8 cm)

Weight: 335 lb (151.9 kg)

Patient Weight Capacity: 500lb (226.7 kg)

5" (12.7 cm) Individual locking swivel casters

MATTRESS

Torso section 2" (5 cm) thick, center and leg sections 3" (8 cm) thick;

Foam rubber pad covered with seamless, durable Naugahyde

Upholstery Color: Graphite

MOTIONS

Height: 23 in to 39 in (58.4 to 99 cm)

Trendelenburg: 0° to ±15°

Fowler Back: 0° to 80°

CONTROLS

Hand Control: Activates height and Trendelenburg motions

Manual Fowler positioning (Optionally motorized)

Foot Control (Optional): Activates height, trendelenburg, and optional motorized Fowler positioning.