

CUBESCAN BIOCON 700



ABOUT

The **CUBEScan BioCon 700** is a Bladder Scanner, which utilizes ultrasound technology to accurately assess bladder volume. The device is portable and user-friendly, allowing for quick and efficient bladder volume measurements using BladderPoint technology without the need for catheterization. It helps healthcare professionals assess urinary retention, monitor urinary output, and manage bladder-related conditions, ultimately improving patient care and comfort.

FEATURES

Quickly and confidently locate the bladder for accurate measurement with BladderPoint technology that puts a bull's eye in the middle of the bladder.

Dependable performance with built-in self-diagnostics to alert you of any hardware or maintenance issues.

Tough polycarbonate casing promotes consistent, uninterrupted service and is supported by a five-year warranty including drop protection.

Easy to use on bariatric patients with rear grip on the probe and no need to fan probe.

DIMENSIONS

Height: 4.1 in (10.3 cm)

Width: 7.95 in (20.2 cm)

Depth: 2.6 in (36.6 cm)

Weight: 2.0 lbs. (0.92 kg)



SPECIFICATIONS

Modes: Male, Female and Child

Storage Capacity: 700 scans

Display: 7" TFT- LCD (WVGA); 16.7M Colors

Range/Accuracy: 0-999 ml; $\pm 15\%$, $\pm 15\text{ml}$

Ultrasound Probe: 3D Sector Scan, 2.6Mhz
Ultrasound Frequency, B-Mode Scan