

MIDMARK RITTER 419

ABOUT

The **Midmark Ritter 419** procedure and exam chair is designed to provide optimal positioning and support for patients during general exams and procedures, particularly those focusing on the shoulder, neck, and head areas. This exam chair has a narrow neck area and a slim backrest, offering a larger work surface for better patient access. The Midmark 419 can support a maximum of 325 pounds and can drop down to a minimum height of 22.4 inches.

FEATURES

Standard table functions: Power height, tilt, back, and foot adjustments with hand or foot control.

Patient positioning in either a chair or flat table position.

Pull-out foot extension.

Height range of 22.4 to 40.4 inches.



SPECIFICATIONS

DIMENSIONS

Height: 22.4 to 40.4 in (56.9 to 102.6 cm)

Width: 23 in (54.4 cm)

Length: 65 in (165.1 cm) (Back down and foot up)

Weight: 515 lb (233.6 kg)

Headrest: Fully extended adds 15 in (38.1 cm)

Footrest Extension: Adjustable up to 8 in (20.3 cm)

MOVEMENTS

Height Range: 22.4 to 40.4 in (56.9 to 102.6 cm)

Foot Section Range: 0° to -90° ± 5°

Table Top Tilt Range: 0° to +30° ± 5°

Back Section Range: 0° to +85° ± 5°

ELECTRICAL / FUSE

Electrical: 115 VAC ±10%, 11 amps, 50/60 HZ

Fuse: 12A, 250V, Time Lag

CLASSIFICATIONS

Class 1, Type B Applied Part, Ordinary Equipment,
Intermittent Operation