

AXIA EXI

ABOUT

The **Axia EXI** Lift combines high performance with sleek design, ensuring safe and effortless patient transfers. This electric-powered, premium transfer system offers a reliable way to move patients securely, minimizing the risk of injury for both patients and caregivers. Tailored for healthcare facilities, the lift is adept at floor transfers and fits beneath most beds in their lowest positions. This makes it an ideal choice for efficiently transferring patients in compact rooms with low beds.

FEATURES

- 6-point spreader bar for various sling options
- Auxiliary up and down lift controls
- 600-pound weight capacity
- Ability to lift a patient from the floor
- 4-inch rear casters with 3-inch front casters



SPECIFICATIONS

DIMENSIONS

Base Length: 47.3 in (120.14 cm)

Base External Width (Opened): 39.3 in (99.82 cm)

Base External Width (Closed): 26.4 in (67.06 cm)

Mast Height: 51.3 in (130.3 cm)

Weight: 126.5 lb (57.38 kg)

MOVEMENTS

Lifting Height (Highest): 68.6 in (174.24 cm)

Lifting Height (Lowest): 21.5 in (54.61 cm)

Base Opened Turning Diameter: 63.4 in (161.04 cm)

Base Closed Turning Diameter: 58.4 in (148.34 cm)

MAX CASTER HEIGHT

Front Caster Height: 4.7 in (11.94 cm)

WEIGHT CAPACITY

Weight Capacity: 600 lb (272.16 kg)