

SCHURE MED EZ LIFT BEACH CHAIR

ABOUT

The E-Z Lift Beach Chair is designed for orthopedic shoulder procedures to ensure the best patient positioning. This chair is lightweight and can accommodate different patient needs, surgeon preferences, and table models.

FEATURES

Greatly decreases the amount of work for staff when moving patients compared to other manual methods.

Articulation range: 0-90°

The sliding headrest adjusts to fit each patient perfectly.

Lateral brace included

Double-ball joint head system safely secures patient's head

Five-piece Deluxe pad set included

Minimum Height: 4 ft 2

Maximum height: 6 ft 4



SPECIFICATIONS

DIMENSIONS

Height: 47 in

Width: 28 in (28.3 cm)

Depth: 30 in (52.7 cm)

Weight: 32 lb (7 kg)

Weight Capacity: 500lbs