

Patient Warmer

STRYKER MISTRAL-AIR MODEL 1100

ABOUT

The **Stryker Mistral-Air 1100** is a patient warmer that works with low-pressure blankets to evenly distribute diffused warm air over the surface of the patient. The Integrated HEPA filter captures 99.99% of particles down to 0.3 microns. The Mistral-Air is a quiet but high-volume blower with over and under-temperature alarms.

FEATURES

Doesn't have individual high-pressure jets of air that can cause the blanket to move.

Integrated HEPA filter.

Default Startup temperature of 100.4°F

3 heat settings.



SPECIFICATIONS

DIMENSIONS

Length: 10.9 in (27.6 cm)

Width: 15.2 in (38.5 cm)

Height: 9.4 in (23.9 cm)

HEAT SETTINGS

Low: 89.6°F (32°C)

Medium: 100.4°F (38°C)

High: 109.4°F (43°C)

GENERAL

Power: 110-120 V 60 Hz 6 amps

Hose: 6 ft (1.8 m)

Power Cord: 13 ft (4 m)