

# STRYKER 1025 SM304 M-SERIES



## ABOUT

The Stryker 1025 SM304 M-Series stretcher combines advanced technology with ergonomic design to enhance patient transport. Featuring Stryker's Zoom motorized drive system, it significantly reduces the physical strain on caregivers, allowing for effortless maneuvering even in tight spaces. This stretcher is built with a robust frame, capable of supporting up to 700 pounds, and includes a power-washable design for easy maintenance. Its versatile features, such as Glideaway siderails and various adjustable positions, ensure patient safety and comfort during transfers. Ideal for high-demand medical environments, the SM304 delivers reliable performance for a range of patient care needs.

## FEATURES

**Zoom Motorized Drive System:** Effortlessly moves the stretcher with minimal manual effort, improving mobility and reducing caregiver strain.

**High Weight Capacity:** Supports up to 700 pounds (317.5 kg) for a wide range of patients, ensuring stability and safety.

**Power-Washable Design:** Built to withstand rigorous cleaning, making it easy to maintain hygiene standards.

**Glideaway Siderails:** Provides a zero-transfer gap to facilitate safe and smooth patient transfers.

**Advanced Positioning Options:** Includes adjustable fowler, knee gatch, and Trendelenburg settings, allowing for versatile patient positioning.



## GENERAL SPECIFICATIONS

**Maximum Weight Capacity:** 700 lbs (317.5 kg)

**Overall Length:** 84 in (213.4 cm)

**Height Range:** 20.75 - 34 in (52.7 - 86.4 cm)

**Fowler Angle:** 0 to 90°

**Knee Gatch Angle:** 0 to 40°

**Trendelenburg / Reverse Trendelenburg:** +16 / -16°

**Mattress:** 3 in (7.6 cm)

**Optional Scale System**

## 26" LITTER MODEL

**Overall Width (Siderails Up):** 33.5 in (85.1 cm)

**Overall Width (Siderails down):** 30.25 in (76.8 cm)

## 30" LITTER MODEL

**Overall Width (Siderails Up):** 37 in (94 cm)

**Overall Width (Siderails Down):** 30.75 in (78.1 cm)