

# STRYKER CORE 2 CONSOLE



## ABOUT

The **Stryker CORE 2 Console** is a versatile and advanced surgical power tool control system. Designed for use in a wide range of procedures, it manages devices like drills, saws, and shavers, providing consistent and reliable performance. With multiple ports for handpieces, footswitches, and irrigation, the console automatically detects connected devices and adjusts settings accordingly. Its user-friendly touch screen interface allows for easy navigation, while the ability to store and export presets simplifies operation across different procedures. The CORE 2 Console's design ensures it meets the demanding needs of surgical environments, enhancing efficiency and safety.



## FEATURES

**Automatic Device Detection:** Instantly recognizes connected devices and adjusts available settings for optimized functionality.

**Touch Screen Interface:** Provides clear visuals and straightforward navigation, enhancing ease of use during procedures.

**Multiple Device Compatibility:** Supports a variety of surgical instruments including drills, saws, and shavers for versatile application.

**Preset Storage and Export:** Allows for the storage and transfer of user preferences, streamlining the setup for recurring procedures.

**Integrated Irrigation Management:** Supports the use of irrigation cassettes, enabling controlled fluid delivery to enhance surgical precision.

## DIMENSIONS

**Height:** 5.4 in (13.72 cm)

**Width:** 13 in (33.02 cm)

**Depth:** 17.4 in (44.2 cm)

**Weight:** 17.3 lb (7.8 kg)

## HIGHLIGHTS

Universal console for use across different specialties.

Variety of compatible instrument portfolios.

Colored illumination rings for matched hand and foot switches.

Can run 2 drills simultaneously.

Future Smart Equipment management compatibility.

# SPECIFICATIONS



## DIMENSIONS

**Height:** 5.4 in (13.72 cm)

**Width:** 13 in (33.02 cm)

**Depth:** 17.4 in (44.2 cm)

**Weight:** 17.3 lb (7.8 kg)

## DISPLAY

**Type:** Touchscreen 24-bit color

**Size:** 7 in (17.78 cm)

**Viewing Angle:** 170°

## PORTS

3 handpiece ports

2 footswitch ports

1 irrigation cassette port

## CUSTOMIZATION

Create and store user profiles

I.D. Touch Software for custom tactile feel and drill response

Adjustable irrigation flow

Adjustable acceleration settings

Detachable electric handswitch option

4 electric foot pedal options

Customizable foot pedal buttons

## HIGHLIGHTS

Universal console for use across different specialties.

Variety of compatible instrument portfolios.

Colored illumination rings for matched hand and foot switches.

Can run 2 drills simultaneously.

Future Smart Equipment management compatibility.