

HILLROM WELCH ALLYN SURETEMP 679



ABOUT

The **Hillrom Welch Allyn SureTemp 679** is a reliable vital signs monitor designed for fast, accurate temperature measurement in clinical settings. It offers advanced thermometer technology for oral, axillary, and rectal readings, ensuring accurate and quick results for patient care. With its intuitive design and large, easy-to-read display, this device helps streamline workflows for medical staff, making it an essential tool for hospitals, clinics, and urgent care settings.

FEATURES

- **Measurement Modes:** Provides oral, axillary, and rectal temperature readings to support a variety of patient needs.
- **LCD Display:** Ensures clear readability with large, easy-to-read digits, even from a distance.
- **Fever Detection Alerts:** Built-in audio alerts notify clinicians of high-temperature readings.
- **Fast Response Time:** Delivers results in as little as 4-6 seconds, speeding up patient care.
- **Ergonomic Handle:** Designed for comfortable use with one hand, reducing operator fatigue.
- **Interchangeable Probes:** Color-coded probes help reduce cross-contamination between patients.
- **Power Management:** Shuts off after a period of inactivity to preserve battery life.



SPECIFICATIONS



DIMENSIONS

Height: 8.1 in (20.6 cm)

Width: 3.6 in (9.1 cm)

Depth: 2.7 in (6.9 cm)

Weight: 1.3 lbs (0.6 kg)

POWER AC

Voltage: 100-240 V AC, 50/60 Hz

Battery Type: 3 AA batteries (for portable use)

ENVIRONMENTAL

Operating Temperature: 50°F to 104°F
(10°C to 40°C)

Storage Temperature: -4°F to 122°F (-20°C to 50°C)

DISPLAY

Screen Type: Monochrome LCD

Backlight: No