

# WELCH ALLYN HILLROM SURETEMP 692 PLUS



## ABOUT

The Hillrom Welch Allyn SureTemp 692 Vital Signs Monitor is an advanced designed to deliver fast and accurate temperature readings in a clinical setting. With its intuitive design and multipurpose functionality, this monitor is ideal for high-volume environments. In addition to temperature measurement, the SureTemp 692 supports oral, axillary, and rectal probes with protective covers to prevent cross-contamination.

## FEATURES

- **Temperature Measurement:** Delivers oral temperature readings in approximately 4 to 6 seconds, reducing patient discomfort and saving time.
- **Multiple Probe Options:** Supports oral, rectal, and axillary probes, with interchangeable covers to minimize cross-contamination.
- **Portable Design:** Compact and lightweight, suitable for mobile and bedside patient monitoring.
- **Basic NIBP Monitoring:** Offers non-invasive blood pressure measurements as an added feature.
- **Patient Monitoring:** Functions as a basic patient monitor for tracking multiple vitals.
- **Monitor Mode:** Provides continuous temperature monitoring for patients who require frequent checks.
- **Customizable Alerts:** Allows configuration of alerts and warnings to notify clinicians of abnormal readings.



# SPECIFICATIONS



## DIMENSIONS

**Height:** 7.0 in (17.8 cm)

**Width:** 2.25 in (5.7 cm)

**Depth:** 3.0 in (7.6 cm)

**Weight:** 10.25 oz (290 g)

## DISPLAY

**Type:** Liquid Crystal Display (LCD)

**Icons:** 3½-digit display with special icons for quick reading

## POWER

**Battery Type:** 3 standard AA batteries

**Battery Life:** Approx. 5,000 readings per set of batteries

## OPERATING ENVIRONMENT

**Operating Temperature:** 60.8°F to 104°F (16°C to 40°C)

**Storage Temperature:** -4°F to 120°F (-20°C to 50°C)

## MEASUREMENT RANGE

**Temperature:** 84.0°F to 108.0°F (28.9°C to 42.2°C)

**Accuracy:** ± 0.2°F in Monitor Mode