

BAXTER HILLROM PROGRESSA



ABOUT

The **Baxter Hillrom Progressa ICU** hospital bed is designed to enhance patient care and comfort while streamlining caregiver workflows. With advanced therapeutic surfaces and continuous lateral rotation therapy, it aids in the prevention and treatment of pressure ulcers and pulmonary complications. The integrated Graphical Caregiver Interface (GCI) provides real-time feedback and controls for easy operation. This bed supports a range of clinical settings, including acute care, critical care, and high-acuity units, and adapts to diverse patient needs with its versatile positions and features. Built with safety and convenience in mind, the Progressa bed promotes efficient care delivery in modern healthcare environments.



FEATURES

Advanced Therapeutic Surfaces: Provides options for continuous lateral rotation and percussion/vibration therapies.

Graphical Caregiver Interface (GCI): Simplifies operation with intuitive controls for bed adjustments, scale, and alarms.

Versatile Bed Positions: Includes chair, Trendelenburg, and FullChair modes for various patient care needs.

Patient Comfort Controls: Includes entertainment, nurse call, and light controls accessible from the bedside.

Enhanced Mobility Options: Equipped with IntelliDrive for powered transport and obstacle detection for safety.

DIMENSIONS

Highest Position (Seat Section to Floor): 35 in (89 cm)

Lowest Position (Seat Section to Floor): 16.5 in (41.9 cm)

Length: 88 in (223.5 cm)

Width: 40.5 in (103 cm)

C-Arm Clearance: 7-15 in (17.8-38.1 cm)

Safe Working Load: 650 lb (295 kg)

Maximum Patient Weight: 500 lb (227 kg)

Radiolucent Sleep Deck (Standard): 17.7 x 23 in (45 x 58.5 cm)

SPECIFICATIONS



DIMENSIONS

Highest Position (Seat Section to Floor): 35 in (89 cm)

Lowest Position (Seat Section to Floor): 16.5 in (41.9 cm)

Length: 88 in (223.5 cm)

Width: 40.5 in (103 cm)

C-Arm Clearance: 7-15 in (17.8-38.1 cm)

Safe Working Load: 650 lb (295 kg)

Maximum Patient Weight: 500 lb (227 kg)

Radiolucent Sleep Deck (Standard): 17.7 x 23 in (45 x 58.5 cm)

MOVEMENTS

Head Section Adjustment: 0° to 67°

Thigh Section Adjustment: 0° to 30°

Foot Section Adjustment: 0° to 70°

Trendelenburg: 13°

Reverse Trendelenburg (Tilt Table Function): -18°

GENERAL

Standard Caster Size: 6 in (15.24 cm)

Cleanability/Infection Control: IPX4 certified

OPTIONAL MATTRESS TYPES

Prevention: Designed to minimize pressure ulcers.

Therapy: Supports advanced therapeutic treatments.

Pulmonary: Optimized for respiratory therapies.

ELECTRICAL SPECIFICATIONS

Voltages Supported: 100, 110, 115, 120, 127, 220, 230, 240 volts