

ARJO SARA COMBILIZER



ABOUT

The **Arjo Sara Combilizer** is a versatile patient positioning and mobilizing aid designed for early mobilization in critical care environments, such as ICUs. Combining the functions of a tilt table, stretcher, and chair, it facilitates rehabilitation by supporting patients in standing, sitting, or supine positions. This multifunctional aid enhances respiratory and cardiovascular function, improves psychological well-being, and reduces the risk of pressure injuries. Built with patient safety and caregiver convenience in mind, the Sara Combilizer is an essential tool for modern critical care and rehabilitation.



FEATURES

Early Mobilization Support: Enables sedated or mechanically ventilated patients to be positioned in standing, sitting, or supine positions for extended periods.

Lateral Tilt Function: Allows the surface to tilt left or right up to 20°, aiding in weight redistribution and rehabilitation exercises.

Enhanced Safety: Includes a manual quick-down function, ensuring patients can be swiftly returned to a supine position in emergencies.

Ergonomic Control Options: Operated via hand controls, a control panel on the handle, and an emergency control box for caregiver ease.

Durable and Mobile Design: Equipped with low-friction casters and brakes, enabling smooth maneuverability in critical care settings.

DIMENSIONS

Length (Seated position): 62.25 in (158 cm)

Maximum Length (Stretcher): 80.5 in (204.5 cm)

Width: 29.5 in (72 cm)

Width of seat: 20.12 in (51 cm)

Depth of Seat: 18 in (46 cm)

Minimum Lifting Height: 23.25 in (58.8 cm)

Maximum Lifting Height: 38.75 in (98.4 cm)

Patient Weight Capacity: 440 lb (200 kg)

Weight: 254 lb (115 kg)

MOVEMENTS

Backwards Tilt: -25° to 0°

Longitudinal Tilt: -25° to +75°

Sideways Tilt (Left and right): 0° to 20°