

MIDMARK 625 BARRIER-FREE EXAM TABLE



ABOUT

The **Midmark Ritter 625** Barrier-Free Exam Table is designed to enhance accessibility, efficiency, and patient care in medical practices. With its low height, advanced positioning options, and integrated technology, it helps reduce physical strain for caregivers while making the exam process more comfortable for patients. The built-in digital scale and connectivity options streamline workflows, allowing for faster and more accurate data collection. Engineered for flexibility, this exam table accommodates a variety of patient needs, including those with limited mobility or bariatric requirements.

FEATURES

Barrier-Free Accessibility: Low minimum height of 18 inches (46 cm) enables easier transfers for wheelchair users and patients with limited mobility.

Advanced Multi-Positioning: Powered seat and back adjustments offer quick and precise patient positioning, including supine, seated, and lithotomy positions.

Premium Comfort: Features elastomeric upholstery for a chair-like feel, enhancing patient comfort during exams.

Ample Storage Options: Equipped with spacious drawers for exam supplies, including optional rear storage for bulky items.

Ergonomic Assist Handles: Provides added stability and security for patients during transfers.

Optional Built-In Digital Scale: Provides accurate weight measurements directly on the table, supporting patients up to 650 lbs (295 kg).



DIMENSIONS AND MOVEMENTS

Height: 18 - 37 in (46 - 94 cm)

Width (Wide-Premium top): 32 in (81 cm)

Width (Ultra-Premium Top): 28 in (71 cm)

Length: 59.5 in (151.13 cm)

Length (With footrest extended): 75.5 in (190 cm)

Weight: 433 lb (197 kg)

Patient Weight Capacity: 650 lb (295 kg)

Paper Roll Holder: Holds 21 x 3.5 in dia. paper roll

Pelvic Tilt: 7°