

STRYKER PRIME X SERIES



ABOUT

The **Stryker Prime X** Series Stretcher is a versatile and innovative platform designed to enhance patient care, imaging accessibility, and mobility in hospital environments. With three available models—1105X, 1115X, and 1125X—the Stryker Prime X provides advanced features tailored to meet diverse clinical needs. Built with Stryker's ClearView Technology, the dual-deck design offers seamless X-ray imaging capabilities while maintaining the durability and mobility required in critical care settings. Whether for transport, imaging, or surgical preparation, the Prime X ensures reliability, safety, and comfort for both caregivers and patients.

FEATURES

ClearView Technology: Dual-deck design allows 360° cassette positioning, delivering high-quality X-ray imaging equivalent to traditional X-ray tables.

Multiple Models for Advanced Mobility:

- **1105X:** Fifth-wheel steering for enhanced maneuverability.
- **1115X:** Big Wheel Advanced Mobility System reduces push force and improves control.
- **1125X:** Zoom Motorized Drive System significantly minimizes manual pushing effort.

Exceptional Weight Capacity: Accommodates patients up to 700 lbs (318 kg), ensuring support for bariatric needs.

Glideaway Siderails: Provide unobstructed lateral imaging access while fully folding for compact storage.

User-Friendly Design: Equipped with dual-sided hydraulic controls, integrated IV poles, oxygen bottle holders, and ergonomic push handles.



MODELS

1105X: Stretcher with Fifth Wheel

1115X: Stretcher with Big Wheel Advanced Mobility

1125X: Stretcher with Zoom Motorized Drive System

DIMENSIONS

Height (1105X and 1115X): 23.25 – 36.5 in (59.1 – 92.7 cm)

Height (1125X): 23 – 34 in (58 – 86 cm)

Width (Siderails up): 38 in (96.5 cm)

Width (Siderails down) (1105X, 1115X, 1125X): 30.5 in (77.5 cm)

Overall Length (1105X and 1115X): 85 in (215.9 cm)

Overall Length (1125X): 86 in (218.4 cm)

Patient Surface: 30 x 75.25 in (76 x 191 cm)

SPECIFICATIONS



MODELS

- 1105X:** Prime X Stretcher with Fifth Wheel
- 1115X:** Prime X Stretcher with Big Wheel Advanced Mobility
- 1125X:** Prime X Stretcher with Zoom Motorized Drive System

DIMENSIONS

- Height (1105X and 1115X):** 23.25 – 36.5 in (59.1 – 92.7 cm)
- Height (1125X):** 23 – 34 in (58 – 86 cm)
- Width (Siderails up):** 38 in (96.5 cm)
- Width (Siderails down) (1105X, 1115X, 1125X):** 30.5 in (77.5 cm)
- Overall Length (1105X and 1115X):** 85 in (215.9 cm)
- Overall Length (1125X):** 86 in (218.4 cm)
- Patient Surface:** 30 x 75.25 in (76 x 191 cm)

POSITIONING

- Backrest:** 0°–90° ($\pm 5^\circ$)
- Trend./Reverse Trend. (1105X and 1115X):** $\pm 17^\circ$ ($\pm 3^\circ$)
- Trend./Reverse Trend. (1125X):** $\pm 16^\circ$

GENERAL

- Weight Capacity:** 700 lbs (318 kg)
- Castor Diameter:** 8 in (20 cm)