

Pain Management Tables


AXIA C-ARM TABLES

ABOUT

The Axia C-Arm Tables surgical table series is intended to be used with mobile x-ray units for minimally invasive procedures on patients as large as 500 lbs.

FEATURES

- Patient Capacity: Axia tables can hold up to 500 lb in patient weight,
- Tabletop Material: Covered with carbon fiber skin/foam core with integral facial cut-out
- Radiolucent Area: 10 square feet of unobstructed Radiolucent area.
- Tabletop: 84"L x 24"W x 1.62"T
- 5" Castors with a full swivel & a wheel lock
- I-Shaped Steel Base with Epoxy Finish
- Patient Straps: Contains 2 patient straps to contain and provide security
- Hand Control: Activates all motions with a click of a button



Axia C-Arm Tables

SPECIFICATIONS

SL1

Movement Modes:

- Height: 30" to 40"

Table Weight: 375 lb (170 kg)

AC Current: 2.5A

SL2

Movement Modes:

- Height: 30" to 40"
- Lateral Tilt: (+ / -) 17°

Table Weight: 390 lb (177 kg)

AC Current: 2.5A

TL3

Movement Modes:

- Height: 30" to 40"
- Lateral Tilt: (+ / -) 17°
- Trendelenburg: (+ / -) 13°

Table Weight: 420 lb (191 kg)

AC Current: 2.5A

STL4

Movement Modes:

- Height: 30" to 40"
- Lateral Tilt: (+ / -) 17°
- Trendelenburg: (+ / -) 13°
- Longitudinal travel: 0" (+/- 5 ")

Table Weight: 590 lb (268 kg)

AC Current: 6.5A

STL5

Movement Modes:

- Height: 30" to 40"
- Lateral Tilt: (+ / -) 17°
- Trendelenburg: (+ / -) 13°
- Longitudinal travel: 0" (+/- 5 ")
- Lateral travel: 10" (+/- 5 ")

Table Weight: 590lbs (268kg)

AC Current: 6.5A

DIMENSIONS

Length: 86 in (218.5 cm)

Width: 26 in (66.0 cm)

Tabletop Length, Width: 84 in (213.3 cm),
24 in (61.0 cm)

I-Shaped Steel Base Length, Width: 82 in (208.3
cm), 26 in (66.0 cm)

Mattress: 2 in (5 cm)

Radiolucent Area Length, Width: 56 in (142.2 cm),
24 in (61.0 cm)

GENERAL

Patient Capacity: 500 lb (227 kg)

Tabletop and Mattress attenuation 1.6mm Al

Patient Straps: 2

Hand Control: Activates all motions

POWER

AC Voltage: 120V~

Frequency: 60 Hz

Electric Shock Protection: Class I

Applied Part Rating: Type B

Mode of Operation: 10% Duty Cycle

1 min on, 9 minutes off