



The **Hillrom Advanta 2 ICU Bed** is designed to help you best meet your daily challenges. Its intuitive and integrated functionality—like positioning, transporting, communication, monitoring, and alerts—deliver the features to help address your top priorities and enhance your experience as a caregiver.

## Features

- IntelliDrive Powered Transport System
- Line-of-Site
- Chair Position
- Care Grip Handle
- Three-Mode Bed Exit System
- Integrated Scale
- One-Step Operation Siderails
- Point-of-Care Siderail Controls
- Electric Vascular Foot



# Specifications

## Physical Specifications

**Frame Width:** 40" (101.6 cm)

### Frame Length with Head / Footboard (w/ Bumpers)

**Normal:** 91" (231 cm)

**Extended:** 94" (238.9 cm)

## Sleep Deck

**Width:** 36" (91.44 cm)

**Normal Length:** 84" (213.4 cm)

**Extended Length:** 87" (221 cm)

## Bed Clearance Under Frame

5.00" (12.70 cm)

## Optional scale system

**Accuracy:** ± 2.2 lbs or 1.0% of patient weight w/e is greater

**Weigh in Any Position:** yes

**Scale Repeatability:** 0.2% of patient weight

## Bed Clearance Under Frame

5.00" (12.70 cm)

## Low Position Floor to Deck

15.75" (42.55 cm)

## High Position Floor to Deck

32.5" (82.55 cm)

## Max Head Elevation

65°

## Max Trend / Reverse Trend

16°

## Caster Size

6" (15 cm)

## Safe Working Load

500 lbs

