



Soma offers the Integra MAYFIELD Triad A1108 Skull Clamp refurbished up to 50% below the cost of new. The MAYFIELD Triad A1108 Skull Clamp offer 3-point skeletal fixation, holding the head and neck in a specific position for Cranial, Spinal or Neurosurgical procedures. The three Skull pins can be adjusted to apply specific clamping force at each point. The Skull clamp offers a sturdy locking mechanism and quick release to ensure the patient's position will not move during the procedure and for quick removal. The skull clamp has the ability to rotate 360 degrees with a swivel rocker arm when force is applied. The skull clamp sits slightly lower than the position of the headrest for better support of the patient.

## Features

- Flexibility of the positioning of the skull pins
- 3-point skeletal fixation
- Offer 360 rotation with a swivel rocker arm
- Quick release lock and easy removal

