



The Oakworks CFPM300 Pain Management C-Arm Imaging Table is a new, cost-effective table designed to accommodate mobile and ceiling suspended C-Arm systems. The CFPM300 is a 3-movement table that has all of the features required for successful pain management, image-guided procedures. Its dynamic height range of 26" to 44" (66cm to 112cm) accommodates low height patient wheelchair transfers to maximum height C-Arm mobility requirements for full lateral imaging.

Features

- 3 Powered Motions
- Motions: Height, Trendelenburg tilt, Lateral Tilt
- Patient Weight Capacity: 500 lbs (227 kg)
- +15/-12° trendelenburg tilt
- Lateral tilt +/-15°
- 22"x84" (56 cm x 213 cm) integrated headrest top

