



The 3502 EZ Slide surgical table permits full body imaging at the click of a button. With 21" of top slide, the patient can be effortlessly moved throughout the imaging window at the simple command of the hand-held control. The low profile base will not get in the way of movable imaging equipment and the low table height can accommodate seated surgeons. With its impressive 700 lb lift capacity, this table is also an ideal solution for bariatric cases.

Features

- 21" horizontal slide
- 700 lb capacity
- 26-44.875" height
- 20° tilt
- 25° trendelenburg
- 20° reverse trendelenburg
- kidney bridge

