



The **SonoScape X5** is a laptop-style Ultrasound. The SonoScape X5 offers ultrasonic imaging for Abdominal, Cardiac, MSK, Reproductive, Vascular, and Small Parts applications. The system offers a portable color Doppler ultrasound machine that can produce crystal clear imaging while still a small compact unit. The Ultrasound can easily and effortlessly move from room to room. The X5 offers a large library of probes and transducers that the system is comparable with, allowing you to provide a comprehensive diagnosis. The ultrasound system offers a 15.6 -inch High-Resolution LCD monitor, a single probe/transducer port, a built-in battery that can quickly charge, and a WiFi DICOM image transfer connection.

Features

- 15.6-inch high-definition monitor
- Built-in battery support for 3 hours of scanning
- WiFi and bluetooth wireless connection
- Multi-beam processing technology
- Spatial compound imaging
- μ -Scan image processing technology
- Pulse inversion harmonic imaging

